

Ground breaking, Life changing™

# Accessories

Stopcocks	117
High-Pressure Injector Lines	118
Extensions / Connectors / Inflation Systems	
Control Syringes	122
Angiographic Needles	

# Accessories Stopcocks

### **1050 PSI STOPCOCKS**

- High pressure capability
- Large inner lumen
- Smooth rotation
- Contents per case: 25

HANDLES OFF/ON	WAY	DESCRIPTION	CATALOGUE
			NUMBER
"OFF"	3-way	with rotator, male	12-002023
"OFF"	3-way	with fixed luer-lock, male	12-002024
"ON"	1-way	with rotator, male	12-001914
"ON"	1-way	with fixed luer-lock, male	12-002022

<sup>\*</sup> check availability with your local Cordis sales representative

### **High-Pressure Injector Lines**

#### **HP INJECTOR LINES 1200 PSI - NON-BRAIDED**

- Non-braided
- Ideal for all power injector applications, offered in a variety of lengths
- Available with or without High-Pressure Rotating Adaptor
- Contents per case: 25

LENGTH	DESCRIPTION	CATALOGUE
(cm)		NUMBER
25	with rotator	12-000813
50	with rotator	12-001771
75	with luer-lock, male	12-000237
75	with rotator	12-001733
120	with luer-lock, male	12-000608
120	with rotator	12-001734

#### **HP INJECTOR LINES 1200 PSI - BRAIDED**

- Braided
- Combines the high pressure capacity with superior flexibility
- Available with or without High-Pressure Rotating Adaptor
- Contents per case: 25

LENGTH	DESCRIPTION	CATALOGUE
(cm)		NUMBER
25	with luer-lock, male	12-001743
50	with rotator	12-001738
75	with luer-lock, male	12-001741
75	with rotator	12-001737
120	with luer-lock, male	12-001740
120	with rotator	12-001736

### **Manifolds 1050 PSI**

• Continuous large internal diameter (2.3 mm diameter) ideal for special procedures requiring small catheters.

DESCRIPTION	CATALOGUE
	NUMBER
2-port, with rotator, right-handed	12-001759
3-port, with rotator, right-handed	12-001765
3-port, with rotator, right-handed	12-001767

# High-Pressure Injector Lines

### **Pressure Monitoring Lines**

• Clear and high durometer for accurate measurements and optimum dynamic response

LENGTH	DESCRIPTION	CATALOGUE
(cm)		NUMBER
15	male/female	12-000235
30	male/female	12-000836
90	male/female	12-001776
120	male /female	12-000236
150	male/female	12-000210

# Extensions / Connectors / Inflation Systems

#### **CATHETER EXTENSIONS**

- Flexible, translucent polyurethane
- Braided
- Designed for 1000 psi
- Luer lock hubs
- Male-female
- 5 units per package



LENGTH	INNER DIAMETER	CATALOGUE
(cm)	(mm)	NUMBER
45	1.8	502-100D
65	1.8	502-101D
100	1.8	502-102D

### **VALVES & HAEMOSTATIC Y-CONNECTORS**

- Available for different wire diameters
- Smooth turning Tuohy-Borst valve
- With rotating adaptor
- Available as a PTCA valve accessory kit with either a tool or a puncture device and a torque device

CASE CONTENT	DESCRIPTION	CATALOGUE NUMBER
25	Large Bore haemostatic Y-connector for catheter up to 9F (Torque device - Guidewire Insertion Tool - Haemostatic Y-connector for Catheters up to 9F)	12-001915
25	Haemostatic Y-connector for catheters up to 9F	12-001855
25	Torque device	12-001624
25	Twin interchange Y-connector for catheters up to 9F	12-001856

#### INFLATION SYSTEMS

INITERITOR STSTEMS	CASE CONTENT	DESCRIPTION	CATALOGUE NUMBER
	6	Disposable Inflation Device 30ml pressure control syringe with manometer measuring pressures between 0 and 30 atm (incl 3-way stopcock)	12-003704

# Extensions / Connectors / Inflation Systems

### **CONTRAST SETS**

CASE CONTEN		CATALOGUE
		NUMBER
	Standard bore tubing	12-002002
	Small bore tubing	12-000303

# **Control Syringes**

### **CONTROL SYRINGES**

- Smooth operation
- Clear, glass-like plastic syringe for clarity
- Contents per case: 25



VOLUME	DESCRIPTION	STOPPER	CATALOGUE
(ml)		(ml)	NUMBER
10	with luer-lock, male	0.5	12-000249
10	with rotator, male	0.5	12-000304

# Angiographic Needles

### **ONE-PIECE NEEDLES**

- Thin wall cannula
- Luer hub
- 10 units per package

CANNULA OD	DESCRIPTION	MAXIMUM GUIDEWIRE (")	LENGTH (cm)	CATALOGUE NUMBER
18G	Thin wall	0.038	7	12-004943
19G	Thin wall	0.032	3.8	12-004925
21G	Thin wall	0.021	3.8	12-004924

### **ONE-PIECE NEEDLES**

- Thin wall cannula
- Luer hub
- 25 units per package

CANNULA OD	DESCRIPTION	MAXIMUM GUIDEWIRE (")	LENGTH (cm)	CATALOGUE NUMBER
18G	Thin wall with snap-on wings	0.038	7.0	502-652
19G	Thin wall with snap-on wings	0.032	7.0	502-654*
20G	Thin wall with snap-on wings	0.021	5.0	502-655*
20G	Thin wall with snap-on wings	0.025	5.0	502-659*

#### **TWO-PIECE NEEDLES**

- Thin wall cannula
- Luer hub

- Solid stylet
- 10 units per package

CANNULA OD	DESCRIPTION	MAXIMUM GUIDEWIRE (")	LENGTH (cm)	CATALOGUE NUMBER
18G	Thin wall	0.038	7.3	12-002131

• 25 units per package

CAN	NULA OD	DESCRIPTION	MAXIMUM GUIDEWIRE (")	LENGTH (cm)	CATALOGUE NUMBER
	18G	Thin wall with snap-on wings	0.038	7.0	502-657
	19G	Thin wall with snap-on wings	0.035	7.0	502-658

Search

# Angiographic Needles

### **250 PSI MANIFOLDS**

		NUMBER
ON	2-port, with rotator, right-handed	12-004948
OFF	2-port, with rotator, right-handed	12-004951
OFF	3-port, with rotator, right-handed	12-004952
	550 PSI Manifolds 3-port	12-001728

### **FLUID ADMINISTRATION SETS**

- Allows for smooth and bubble-free of closed system
- Vented or non-vented
- Inner diameter 3 mm

DESCRIPTION	CATALOGUE	
	NUMBER	
standard tubing, vented with drip chamber	12-000214	

DESCRIPTION

Search

CATALOGUE